FOUNDATIONONE® CDx

REPORT  
EU

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
**FOUNDATIONONE® CDx**

**ABOUT THE TEST** FoundationOne®CDx is a next-generation sequencing (NGS) based assay that identifies genomic findings within hundreds of cancer-related genes.

**PATIENT**  
 DISEASE Lung adenocarcinoma  
 NAME Not Given  
 DATE OF BIRTH Not Given  
 SEX Female  
 MEDICAL RECORD # Not Given

**PHYSICIAN**  
 ORDERING PHYSICIAN Not Given  
 MEDICAL FACILITY Not Given  
 ADDITIONAL RECIPIENT Not Given  
 MEDICAL FACILITY ID Not Given  
 PATHOLOGIST Not Given

**SPECIMEN**  
 SPECIMEN SITE Not Given  
 SPECIMEN ID Not Given  
 SPECIMEN TYPE Not Given  
 DATE OF COLLECTION Not Given  
 SPECIMEN RECEIVED Not Given

**PATIENT Sample** Jane  
**TUMOR TYPE** Lung adenocarcinoma  
**REPORT DATE** 01 Jan 2018  
**GRF#** XXXXXXXX

**Biomarker Findings**  
 Tumor Mutational Burden - TMB-Intermediate (11 Muts/Mb)  
 Micros

**Genomic Findings**  
 For a complete list of the genes assayed, please refer to the Appendix.  
**EGFR** amplification, L858R  
**PTCH1** T416S  
**CDKN2A/B** loss  
**RBM10** Q494\*  
**TP53** R267P

**7 Disease relevant genes with no reportable alterations: KRAS, ALK, BRAF, MET, RET, ERBB2, ROST**

**14** Therapies with Clinical Benefit  
 0 Therapies with Lack of Response

**18** Clinical Trials

BIOMARKER FINDINGS	THERAPIES APPROVED IN THE EU (IN PATIENT'S TUMOR TYPE)	THERAPIES APPROVED IN THE EU (IN OTHER TUMOR TYPE)
<b>Tumor Mutational Burden - TMB-Intermediate (11 Muts/Mb)</b> 9 Trials see p. 14	Atezolizumab Durvalumab Nivolumab Pembrolizumab	Avelumab
<b>Microsatellite status: US - MS-Stable</b>	No therapies or clinical trials. see Biomarker Findings section	
<b>GENOMIC FINDINGS</b> <b>EGFR - amplification, L858R</b> 4 Trials see p. 16	Afatinib Erlotinib	Cetuximab Lapatinib
<b>PTCH1 - T416S</b> 5 Trials see p. 17	Osimertinib None	Sonidegib Vismodegib

**GENOMIC FINDINGS WITH NO REPORTABLE THERAPEUTIC OR CLINICAL TRIALS OPTIONS**

For more information regarding biological and clinical significance, including prognostic, diagnostic, germline, and potential chemosensitivity implications, see the Genomic Findings section.

**CDKN2A/B - loss** ..... p. 5    **TP53 - R267P** ..... p. 6  
**RBM10 - Q494\*** ..... p. 5

NOTE: Genomic alterations detected may be associated with activity of certain approved therapies; however, the agents listed in this report may have varied clinical evidence in the patient's tumor type. Therapies and the clinical trials listed in this report may not be complete and exhaustive. Neither the therapeutic agents nor the trials identified are ranked in order of potential or predicted efficacy for this patient, nor are they ranked in order of level of evidence for this patient's tumor type. This report should be regarded and used as a supplementary source of information and not as the single basis for the making of a therapy decision. All treatment decisions remain the full and final responsibility of the treating physician and physicians should refer to approved prescribing information for all therapies.

Therapies contained in this report may have been approved by the US FDA.

Electronically Signed by Julia A. Elvin, M.D., Ph.D. • Jeffrey S. Ross, M.D., Medical Director • 30 November 2017  
 Foundation Medicine, Inc. • 1-888-988-3639

Sample Preparation: 150 Second St., 1st Floor, Cambridge, MA 02141 • CLIA: 22D2027531  
 Sample Analysis: 150 Second St., 1st Floor, Cambridge, MA 02141 • CLIA: 22D2027531

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10. Schrock AB et al. Clin Cancer Res 2016; 22:3281-3285.
11. Rozenblum AB et al. J Thorac Oncol 2017; 2:258-268.